



以下為常規在庫品，均可非標準，均可依樣定制
 Below are regular stock products, we can also make non-standard products according to customers' drawings or samples



部分尺寸 some size

直徑 ΦD Diameter	2.0 (凹頭)	2.5 (凹頭)	3.0 (凹頭)	3.5 (凹頭)	4.0 (凹頭)	2.0 (圓頭)	2.5 (圓頭)	3.0 (圓頭)	3.5 (圓頭)	4.0 (圓頭)
柄徑 ΦD3 Shank Dia.	8.5	8.5	8.5	8.5	8.5	8.5	8.5	8.5	8.5	8.5
總長 L Overall length	168	168	168	168	168	168	168	168	168	168



限位套材料:
 可不同顏色
 高耐磨醫
 塑料類
 可不同顏色
 高品級
 屬類

Material of stop collars:

- ① Can be high wear-resistant medical plastic with different color.
- ② Can be high grade metal with different colors.

