



醫用工具及配件系列 (鑽針、骨鑽、骨科工具系列)
 Medical tools and parts series (Drill points,
 Orthopedic drills, Orthopedic tools series)
 医療器具及び部品シリーズ (ドリルピン、骨ドリル、整形外科工具)


(超防銹型(AA)&超高耐磨型(A))
 (High anti-rust type AA)
 & super wear-resistant type A)
 (超防銹型(AA)&超耐磨型(A))




圓型多刃
Circular multi-flute
丸型多刃
Ball-flute
ボール型刃




任意深淺的容屑槽及刃數、螺旋角..... 270°
CNC成型研磨。
It can be any depth of chip flute with any numbers of flute and helix angle..... 270° profile grinding.





成型錐度鉸型 Type of profile taper reamer
成型コーニング絞り型






鑽鉸型 Type of drilling-reamer
ドリルリーマ型







A向 Direction A
A方向

可超硬研磨任意螺旋角、刃口的鋒利程度類型要求
We can super hard grind any helix angle and satisfy your requirement of any degree of sharpness for cutting edge.






成型、錐度鉸型 Profile, taper-milling
成型、コーニング・フライス型



★ 依圖依樣尺寸/公差/要求復成型加工，可生產極為複雜、精密的柄部形狀。Compound contour machining according to the sizes/tolerances/requirements of drawings and samples, we can produce extreme complex and accurate holder profile.
 圖面、サンプル、寸法または公差の要求により復成型加工を行い、極めて複雑、精密な柄部の形状を生産できる。



—243—