

醫用工具及配件系列 (鑽套、襯套)
 Medical tools and parts series (Drill stopper, sleeve)
 医療器具及び部品シリーズ (ドリルヘッドストッパー、スリーブ)

(可不同顏色處理&高精密成型)
 (We can make different colors & high precision shaping)
 (多色處理&高精密成型可)





定位夾具 (可以加工任意規格尺寸及鬆緊度
 (不銹鋼/鈦/鈦合金)
 Positioning fixture (We can process any
 size and tightness (stainless steel
 /titanium/titanium alloy)



- 定位夾具方式:
- ①、卡口型 (方便、傳統型一般的卡緊程度)
 - ②、螺紋型 (新型、高轉速不易脫落, 但鑽頭也要配合螺紋)

- Positioning fixture type:
- 1. Bayonet type (convenient and traditional type, general chucking level)
 - 2. Thread type (new type, not easy to fall off at high rotary speed, but drills need to match with thread)



★ 超高防銹高硬度高耐磨、高剛性、高抗衝擊、高韌性不銹鋼、鈦、鈦合金等高精度、超細、超長、超硬加工成型。
 Ultra-high rust-resistant, high-hardness, high-wear-resistant, high-rigidity, high-impact-resistant & high-toughness stainless steel, titanium, titanium alloy, etc., with high-precision, ultra-fine, ultra-long, super-hard machine shaping.