

醫用工具及配件系列 (手術刀/克氏針/精微切刀)

Medical tools and parts series (Surgical knives /
Kirschner wires / Subtle cutter)

醫療器具及び部品シリーズ(カッター、メス/キルシュナー鋼線/マイクロカッター)

(超鋒利刃口& 超高耐磨& 超高防銹& 高精密度)
(Ultra-sharp edge & ultra-high wear-resistant
& ultra-high rust & high precision)
(超銳利刃先&超耐磨耗&超防銹&高精度)

CNC磨削中心,可超精研磨任意形狀,結構,精度公差的超鋒利
成型刃口.超高防銹不銹鋼(AA)(HRC54° ± 2° 斷面)

CNC Grinding Center, which can superfinish grind the ultra-sharp profile edge
with any shape, structure and precision tolerance, also the ultra-high anti-rust
stainless steel (AA) (HRC54° ± 2° Section)

精研磨超鋒利刃
(非機加工成型的)
Precision grinding super sharp
(non-machining molding)

精研磨,超鋒利刃(非機加工成型的)
Precision grinding super sharp
(non-machining molding)

全精研磨之平斜面
Flat bevel of full
finish grinding

全精研磨之微
凹斜面,刃超鋒利
Slightly concave bevel
of full finish grinding,
ultra-sharp edge

可微細、成型、特殊形狀、
超鋒利刃CBNC、
高光潔度硬研磨成型

It can be micro, forming,
special-shaped, super
sharp-edged CBNC
and high finish hard
ground to mold