

醫用工具及配件系列 (手術刀/克氏針/精微切刀)

Medical tools and parts series (Surgical knives / Kirschner wires / Subtle cutter)


醫療器具及び部品シリーズ(カッター、メス/キルシュナー鋼線/マイクロカッター) (超銳利刃先&超耐磨耗&超防鏽&高精度)

(超銳利刃口&超高耐磨&超高防鏽&高精度)  
(Ultra-sharp edge & ultra-high wear-resistant & ultra-high rust & high precision)

依圖依樣尺寸公差要求可生產特殊、複合、異型、超大超細超長組合成型的系列刀具。

We can produce special, compound, special-shaped, super-large, super-thin/super-long combined profile series tooling according to requirements of drawings, samples and dimension tolerance.

圖面、サンプル、または寸法の公差要求により特殊、複合、異型または超大、超細、超長の組合せ成型のシリーズカッタを生産できる。



可加工任意的槽型及深度、精度要求

We can process with any groove type and any depth, accuracy requirement

★可超硬、超精研磨。任意的柄部類型。芯厚。容屑槽寬度。刃背方式。鑽尖方式。內冷卻孔方式。刃部結構。刃口角度。刃口鋒利程度之要求。及高可至±0.003mm (3µm) 的寸公差要求 We can produce with super hard and super finish grinding and satisfy your requirement of any type of shank, web thickness, width of chip flute, land, drill tip, inner cooling hole, structure of cutting edge, angle of cutting edge, degree of sharpness of cutting edge, and the dimension tolerance can be up to ±0.003mm (3µm)

