


材料:

- ①超高防銹高耐磨不銹鋼系列(AA)(HRC54° ± 2° 斷面)(常備常規)
- ②一般防銹高耐磨不銹鋼系列(A)(HRC45° ~64° 斷面系列)



Material:

- ①Super-high anti-rust high wear resistance stainless steel series(AA)(HRC54° ± 2° section)(regularly stock)
- ②General anti-rust high wear resistance stainless steel series(A)(HRC45° ~64° section series)

部分尺寸(AA):some size (AA)

	刃徑 $\phi d1$ Diameter	M3.75	M4.05	M4.25	M4.55
	刃長 $L1$ Length of cut 刃の長さ	8	8	8	8
	柄徑 $\phi d$ Shank Dia.	2.35	2.35	2.35	2.35
	總長 $L$ Overall length 總長さ	25.5	25.5	25.5	25.5

- ★、可研磨任意規格，錐度、刃徑;刃長;總長;定位;色圈……及尺寸公差要求。
- ★、We can grind any specification, taper;flute diameter;flute length; total length;location,color code and size tolerance as per the requirement.



部分尺寸(AA):some size (AA)

	刃徑 $\phi D1$ Diameter	M3.5	M4.1	M4.7
	刃長 $L1$ Length of cut 刃の長さ	6.5	6.5	6.5
	柄徑 $\phi d$ Shank Dia.	2.35	2.35	2.35
	總長 $L$ Overall length 總長さ	32.5	32.5	32.5

\* 必要可帶穿透&不穿透，不同等級加工的中空、空芯內孔加工。

With or without penetration, hollow & hollow core inner hole process with different grades.  
 必要に応じて、貫通付き&貫通なし、各レベル加工での中空、中空内穴加工。