




醫療工具及配件系列 (牙科銑銼/氧化鋯車針)
 Medical tools and parts series (Dental milling bur/Zirconia turning pin)
 醫療器具及び部品シリーズ (齒科フライス削り・やすり/酸化ジルコニウムバー)

超高耐磨型&高抗衝擊型
 Super high wear resistance & high impact resistance
 超耐磨耗型&高耐衝擊型

刃徑 φD Diameter	刃長L1 Length of cut 刃の長さ	有效長L2 Effective length	柄徑 φ d Shank Dia.	總長L Overall length 總長さ
0.6	3.5	4.5	4.0	47
0.6	4.0	15	3.175	47/37.5
0.6	4.0	12	3.0	47/37.5
1.0	3.0	17	4.0	47
1.0	4.4	14.5	3.0	47/37.5
1.0	4.0	16	4.0	47
1.0	4.0	12	3.0	47/37.5
1.0	5.0	15	3.0	57/37.5
1.0	5.0	17	3.0	47/37.5
1.0	6.6	17	3.0	37.5
1.0	4.0	17	3.175	47/37.5
1.0	6.0	15	4.0	47/37.5
2.0	10	17.5	3.0	35
2.0	5.0	16	3.0	47/37.5
2.0	6.0	15	3.0	47/37.5
0.6	1.6	9.0	3.0	35
1.0	2.0	16	3.0	35
2.0	10	16	3.0	35



刃徑 φD Diameter	刃長L1 Length of cut 刃の長さ	有效長L2 Effective length	柄徑 φ d Shank Dia.	總長L Overall length 總長さ
2.0	6.0	20	3.175	47/37.5
2.0	6.0	15	4.0	47
2.0	8.0	15	4.0	47
2.0	6.0	15	4.0	47
2.0	10	18	3.0	57/37.5
2.0	4.0	16	4.0	47
2.0	6.0	17	4.0	47
2.0	6.7	17.6	3.0	47/37.5
2.0	10	18	3.0	57/37.5
2.0	6.0	8.0	3.0	47
2.0	11.5	20	3.175	47/37.5
3.0	8.0	7.0	4.0	47
3.0	8.0	15	4.0	47
3.0	8.0	12	3.0	47/37.5
~	~	~	~	~

若需其它尺寸，可以依要求生產！ If you need other sizes, we can also customize for you according to your demands!

※常規總長為47.00mm及37.50mm

※Generally total length is 47.00mm or 37.50mm

可以依圖依樣鏡面研磨加工任意材種，形狀、特殊、異型、精度公差之產品。 We can make tools with any material, any shape, special kind, special type, abnormal shape and precision and tolerance by mirror surface grinding according to your drawings or samples