

醫用工具及配件系列(內冷却孔鑽頭)

Medical tools and parts series(Drill with internal cooling hole)

医療器具及び部品シリーズ (内冷却穴ドリルヘッド)

(超耐腐蝕高耐磨(AA)&高精特殊非標成型)

(Ultra-high anti-rust(AA)&high precision special nonstandard forming)

(超耐摩擦型/高耐衝擊用/高精度)

★、各種鑽頭鑽套可以採用任意材料及幾何形狀、成型超精、超硬研磨

Every drill and stopper can be made with any material and to be ground to any geometrical shape by super finish and super hard grinding

※ 必要可帶穿透&不穿透，不同等級加工的中空、空芯內孔加工。

With or without penetration, hollow & hollow core inner hole process with different grades.

必要に応じて、貫通付き&貫通なし、各レベル加工での中空、中空内穴加工。

种植牙顺序图  
Implant sequence diagram



★可超硬、超精研磨。任意的柄部類型、芯厚、容屑槽寬度、刃背方式、鑽尖方式、內冷却孔方式、刃部結構、刃口角度、刃口鋒利程度之要求，及高可至±0.003mm (3μm) 的尺寸公差要求 We can produce with super hard and super finish grinding and satisfy your requirement of any type of shank, web thickness, width of chip flute, land, drill tip, inner cooling hole, structure of cutting edge, angle of cutting edge, degree of sharpness of cutting edge, and the dimension tolerance can be up to ±0.003mm (3μm)