


醫用工具及配件系列 (扭力扳手/转接头)
 Medical tools and parts series (Torque Wrench/Adapter)
 醫療器具及び部品シリーズ (トルクレンチ/アダプター)

(可不同顏色處理&複雜成型加工)
 (different colors treatment & complex forming processing available)
 (多色處理&複雜成型加工可)



可配置不同顏色及
不同材料的圓圈帽

We can do round caps with
different colors and materials

可直接配用CFDA
標準手機柄
Could directly match
with the CFDA
handpiece shank

可配鋼圈卡更緊式
Can be fitted with
steel ring (chucking
tighter type)

依要求可定制超小
或超大扭力型
We can mark the scale
value as customers
require

閉口型
(柄帶內冷卻、通孔)
Closed type (shank with inner
cooling and through-hole)

可依要求標刻度值
We can mark the scale
value as customers require

凹配concave matching

凸配convex matching

開口型
(柄不帶內冷空)
Open type (shank
without inner cooling)

可黃色/黑色/
光亮/亞光
之處理
Yellow/Black
/Shiny/Matt
treatment

★ 超高防銹高硬度高耐磨、高剛性、高抗衝擊、高韌性不銹鋼、鈦、鈦合金等高精度、超細、超長、超硬加工成型。
 Ultra-high rust-resistant, high-hardness, high-wear-resistant, high-rigidity, high-impact-resistant & high-toughness stainless steel,
 titanium, titanium alloy, etc., with high-precision, ultra-fine, ultra-long, super-hard machine shaping.