

醫用工具及配件系列（手術刀/克氏針/精微切刀）

Medical tools and parts series (Surgical knives / Kirschner wires / Subtle cutter)

医療器具及び部品シリーズ(カッター、メス/キルシュナー鋼線/マイクロカッター)

(超鋒利刀口 & 超高耐磨 & 超高防銹 & 高精度)
(Ultra-sharp edge & ultra-high wear-resistant & ultra-high rust & high precision)
(超銳利刃先 & 超耐磨耗 & 超防銹 & 高精度)

CNC磨削中心,可超精研磨任意形狀,結構,精度公差,的超鋒利成型刀口.超高防銹不銹鋼 (AA) (HRC54° ± 2° 斷面)



CNC Grinding Center, which can superfinish grind the ultra-sharp profile edge with any shape, structure and precision tolerance, also the ultra-high anti-rust stainless steel (AA) (HRC54° ± 2° Section)

全精研磨,超鋒利刃(非機加工成型的); 可微凹斜面, 刃超鋒利
Full finishing, ultra-sharp edge (not machined);
The bevel can be slightly concave, the edge is super sharp

可微細、成型、特殊形狀、超鋒利刃CBNC、高光潔度硬研磨成型

It can be micro, forming, special-shaped, super sharp-edged CBNC and high finish hard ground to mold

超鋒利全精研磨刃口
Super sharp full-finish grinding edge



批量可協助開模注塑.
If in batch, we can assist in mold injection

