


不同形狀之頭部 (部分祇供參考用, 實際可依樣依圖形狀, 尺寸公差生產)
 Heads of various shapes (Just for reference, actual production as per the samples,
 pictures and size tolerance)



可以生產不同的頭部形式、形狀、精研磨成型。
 We can produce various head types, shapes and
 fine grind to mold.

可依要求生產不同頭部材料類型 (陶瓷類、鎢鋼類、不銹鋼類、
 鈦及鈦合金類……等系列)
 We can produce heads from various types of material as per the
 requirement, like ceramic, carbide, stainless, titanium and
 titanium alloy, etc.



刃部材料類型:

1. 超高防銹高耐磨高硬度不銹鋼類(AA)(HRC54° ± 2° 斷面).(常規常備)
2. 常規標準防銹高耐磨不銹鋼類(A)(HRC45° ~64° 斷面)系列.
3. 鈦及鈦合金類(HRC35° ~42°)系列.
4. 陶瓷類 (HRA85° ~94°)系列, (HRC67° ~84°)系列.
5. 鎢鋼類 (HRA81° ~94°)系列, (HRC60° ~84°)系列.

Cutting edge material:

1. Ultra high rust-proof, high wear-resistant & high-hardness stainless steel series (AA)(HRC54° ± 2° Section).(regularly stock)
2. General anti-rust & high wear-resisting type (A) (HRC45° ~64° Section) series.
3. Ti and Ti alloy (HRC35° ~42°) series.
4. Ceramic (HRA85° ~94°) series, (HRC67° ~84°) series.
5. Carbide (HRA81° ~94°) series, (HRC60° ~84°) series.